

The Serine Protease Matriptase 2 Inhibits Hepcidin Activation by Cleaving Membrane Hemojuvelin

Laura Silvestri, Alessia Pagani, Antonella Nai, Ivana De Domenico, Jerry Kaplan, and Clara Camaschella

Supplemental Experimental Procedures

Immunofluorescence analysis

HeLa cells were grown for 16 hrs before transfection on glass cover slips in 12-well plates. Transfections were carried out using 1 µg (microgram) of plasmid DNA and 3 µl of the liposomal transfection reagent Lipofectamine in 1 ml of OptiMem, according to the manufacturer's instructions. After 12 hrs, the medium was replaced and, 24 hrs later, cells were washed with PBS and fixed in 4% paraformaldehyde/PBS for 30 min at room temperature, followed by three PBS washes. Cells were permeabilized by incubation in PBS containing 0.2% Triton X-100 and 10% donkey serum for 1 hr. The cells were then incubated for 2 hrs at room temperature with rabbit polyclonal anti-FLAG (1:250 dilution; Santa Cruz). After PBS washing, primary antibodies were visualized using the donkey anti-rabbit conjugated with Alexa Fluor 594 (Molecular Probes, Eugene, OR). Cover slips were mounted on glass slides using Dako Mounting Medium (Dako, Carpinteria, CA) and visualized at a 63X magnification under confocal microscope BioRad Leica TCS SP2 AOBS (Leica Microsystems, Bannockburn, IL).

Cloning strategies

Mutagenesis on matriptase-2 wild type was performed using the following oligonucleotides:

R774C sense: 5' GGCCTCAGTGGCTGGTTCCTGGC 3'

R774C antisense: 5' GCCAGGAACCAGCAGCCACTGAGTGCC 3'

The matriptase-2^{MASK} variant, lacking the serine protease domain (Du et al, 2008), was obtained by amplification from the wild type cDNA by using the following primers:

MASK-sense: 5' GGAATTCTGAGCTGTCTCGGC 3' (EcoRI tail)

MASK-antisense:

5'GCTCTAGATCACTTGTTCATCGTCGTCCTTGTTAGTCAAGCTTGCCACAGTC
ACAGTGC 3' (FLAG, XbaI tail).

The amplicon was cloned into the EcoRI and XbaI sites of pcDNA3.1(+).

Expressing vector encoding uromodulin cDNA (UMOD, cMYC tagged) was a gift of Dr. Luca Rampoldi (Dulbecco Telethon Institute, Molecular Genetics of Renal Disorders, San Raffaele Scientific Institute, Milan, Italy).

Incubation of s-HJV with matriptase-2-expressing cells

HeLa cells were transiently transfected for 18 hrs with HJV^{wt} expressing vector as described above, and then incubated in serum-free media. After 24 hrs, media were collected, centrifugated for 10 min at 800 g to eliminate detached cells and filtered at 0.22 µm.

Matriptase-2^{wt} (MT2^{wt}) or mock- transfected cells were then incubated with the HJV-derived culture media. After 24 hrs, media were collected, concentrated using 5 kDa molecular weight (MW) cutoff ultrafiltration (Amicon Ultra; Millipore, Billerica, MA), and analyzed by SDS-PAGE.

Incubation of s-matriptase-2 with HJV-expressing cells

HeLa cells were transiently transfected for 18 hrs with matriptase-2^{wt} (MT2^{wt}) or empty vectors, and then incubated in serum-free media. After 24 hrs, media were collected, centrifugated for 10 min at 800 g to eliminate detached cells and filtered at 0.22 μm .

HJV^{wt} transfected cells were then incubated with the MT2- or mock-derived culture media. After 24 hrs, media were collected and concentrated using 5 kDa molecular weight (MW) cutoff ultrafiltration (Amicon Ultra; Millipore, Billerica, MA). Cells were incubated with PI-PLC as described above for m-HJV analysis, and then lyzed. Whole cellular extracts, PI-PLC-derived samples and media were then analyzed by SDS-PAGE.

Legend to Figures

Figure S1. Immunofluorescence analysis of wild type and mutants matriptase-2. HeLa cells were transiently transfected and analysed for immunofluorescence. Matriptase-2 (MT2) was revealed by anti-FLAG antibody. UP: unpermeabilized cells; P: Triton X-100 permeabilized cells. The white Bar represents 10 μm .

Figure S2. HeLa cells were cotransfected with the GPI-anchored protein uromodulin cDNA and matriptase-2 expressing plasmid, or the empty vector. Whole cellular extracts, concentrated media and PI-PLC-derived supernatants were analyzed by Western blot. Anti-cMYC was used to detect uromodulin. s-UMOD indicates soluble and m-UMOD membrane-associated uromodulin. The equal loading was verified by α -tubulin.

Figure S1




Figure S2

