



Supplemental Fig 2



Supplemental Figure 2. Effect of PTK6 kinase inhibitor 4f (Millipore Sigma) on SNAIL levels in TNBC cells

A) Validation of 4f as a PTK6 kinase inhibitor. MCF10A cells overexpressing a constitutively active PTK6 were treated with the indicated concentrations of 4f. Lysates were probed with anti-phosphoPTK6 (tyrosine 342) antibody. B) MDA-MB231 cells were treated with 4f and levels of SNAIL were assessed in lysates.