

Table S2.**A**

fibril (μg)	<i>N</i> (Total)
40	5
10	5
4	4
2	4
1	4
0.1	3
0.01 (10 ng)	3
0.001 (1 ng)	4
0.0001 (0.1 ng)	4
0.00001 (10 pg)	4

B

Treatment	<i>N</i> (Total)
monomer-saline	3
saline	4
100°C 3min	4
120°C 20min	4
120°C 20min/0.1% SDS	3
134°C 20min	3
134°C 20min/0.1% SDS	3

C

brain samples	<i>N</i> (Total)	3M	9M
control	4	2	2
DLB	4	2	2
MSA-2 Pu	5	5	not tested
MSA-2 Pu 100°C3min	5	5	not tested
MSA-2 Pu 134°C20min	4	4	not tested